

Introduction

Bispecific antibodies represent a transformative approach in next-generation biologics, offering enhanced therapeutic efficacy through dual-targeting capabilities. Despite their promise, the efficient and rapid identification of optimal lead candidates remains a significant challenge in drug discovery. The Invenra B-Body® Platform addresses these obstacles by providing an integrated, end-to-end solution that accelerates the discovery process. This platform combines rapid de novo antibody identification, seamless bispecific production, and comprehensive functional screening within a streamlined workflow.

Central to workflow success are our proprietary antibody libraries, engineered for exceptional human-like characteristics and developability. These libraries facilitate the high-throughput generation of hundreds of bispecific candidates across a variety of formats, including 1x1, 2x1, and 2x2 configurations. Notably, the B-Body Platform consistently delivers assay-ready, high-quality bispecific molecules in under four months. Moreover, these molecules exhibit developability and stability metrics comparable to monoclonal antibodies, minimizing downstream development risks.


By integrating data-driven diversity and robust expression technologies, the B-Body Platform sets new standards for speed and reliability in bispecific antibody discovery. Here, we describe capabilities to streamline lead selection, reduce development timelines, and enhance the probability of clinical success, positioning the platform as a cornerstone of next-generation biotherapeutic development.

Phage Library Design


Structured to Enhance Natural Diversity

- Comprehensive Germline Representation:** Libraries encompass a diverse array of human germline sequences, carefully chosen to maximize clinical relevance and therapeutic potential.
- Preservation of Native Pairings:** Natural complementarity-determining region (CDR) pairings are maintained, ensuring minimal liabilities and enhanced developability.
- Optimized Germline Frameworks:** Four IGHV germline frameworks are utilized, each selected for their superior biophysical behavior, clinical success rates, and diverse representation across the antibody repertoire.
- Natural Diversification:** Libraries are engineered to retain native CDR-framework pairings, avoiding artificial shuffling and preserving the integrity of natural antibody structure.
- Extensive Diversity Coverage:** Over 30 specialized sub-libraries expand discoverable diversity, providing a robust platform for identifying high-affinity, developable leads.

Example Sublibrary



Example Selection Campaign



VH Germlines




VL Germlines



*Residues colored by biophysical property

Natural Diversity of CDRH3 Length Across Sub-Libraries



Humanness Scoring of Library Derived New Binders vs. Benchmarks




Figure 1. LSTM-Based Model for Evaluating Binder Humanness

- Correlation with Clinical Outcomes:** Analysis of clinical antibody datasets reveals a strong correlation between humanness scores and reduced anti-drug antibody (ADA) rates, highlighting the importance of humanness in therapeutic antibody development.
- High Humanness Across Binders:** Both benchmark and newly identified binders surpass the established humanness threshold, with many new binders achieving scores indistinguishable from fully human antibodies (score = 1).

Invenra's Rapid Bispecific Antibody Discovery: Antigens to Novel Lead Bispecifics in 4 Months




Figure 2. Rapid Bispecific Antibody Discovery Workflow

This schematic outlines a streamlined four-step process enabling the discovery of novel bispecific antibody leads within just four months. The workflow is designed to produce bispecific candidates rigorously qualified for affinity, cross-reactivity, developability, and yield. A comprehensive matrix analysis evaluates all possible permutations of lead binder configurations in parallel, ensuring the identification of the optimal Fab arm arrangement for enhanced therapeutic potential.

The B-Body Platform: Enabling Rapid Bispecific Antibody Discovery




Figure 3. Overview of the B-Body Platform for Accelerated Bispecific Discovery

The B-Body Platform streamlines bispecific antibody discovery and development by addressing common challenges encountered in the production pipeline. Candidate bispecific antibodies are expressed at high yields and purified efficiently in a single step, enabling rapid evaluation. The platform facilitates a seamless transition from early-stage discovery to cell line development and high-yield CMC (chemistry, manufacturing, and controls), ensuring continuity and scalability within a unified system.

Validation of the B-Body Bispecific Platform Matrix Approach to Antibody Discovery & Optimization

Name	hVGene	lVGene	HEK Titer (mg/L)	DSF (°C)	SDS-PAGE AS100 (mM) J501 (mM)	HIC Retention Time (Min)	SMAC Retention Time (Min)	Poly-Specific, Reag. SMC Score (0-1)	AC-SIMS Δmax (nm) Average
Adalimumab (ADA)	IGHV3-9*01	IGKV1-27*01	134.9	71.0	900.0	8.8	8.7	0.00	1.1
Atezolizumab (ATE)	IGHV3-23*04	IGKV1-NL1*01	164.1	73.5	300.0	13.4	19.3	0.07	15.0
Canakinumab (CAN)	IGHV3-33*01	IGKV6D-21*02	45.7	72.0	800.0	9.3	8.7	0.00	0.7
Dactlizumab (DAC)	IGHV1-46*01	IGKV1-5*01	245.1	74.0	900.0	9.3	8.8	0.00	-0.1
Elotuzumab (ELO)	IGHV3-7*05	IGKV1-27*01	213.2	83.5	700.0	10.3	9.3	0.00	-0.2
Evolocumab (EVO)	IGHV1-18*01	IGLV2-14*01	260.7	65.0	700.0	10.4	9.1	0.20	2.2
Farletuzumab (FAR)	IGHV3-30*03	IGKV1D-33*01	220.8	75.5	800.0	9.5	9.1	0.00	-0.5
Golimimumab (GOL)	IGHV3-30*01	IGKV3-11*01	163.2	70.0	0.0	11.4	12.7	0.23	23.0
Guselkumab (GUS)	IGHV5-101*04	IGLV1-40*01	167.3	69.5	700.0	11.4	9.2	0.47	3.4
Trastuzumab (HER)	IGHV3-11*05	IGKV1-5*05	159.5	78.5	800.0	9.7	8.8	0.00	2.0
Ipilimumab (IPI)	IGHV3-30*01	IGKV3-20*01	169.6	73.0	400.0	11.6	13.0	0.23	10.4
Ixekizumab (IXE)	IGHV1-46*01	IGKV2D-29*02	97.3	83.0	500.0	10.9	9.1	0.81	20.0
Mepolizumab (MEP)	IGHV2-70*20	IGKV4-1*02	221.5	78.5	900.0	9.2	8.8	0.00	-1.0
Ramucirumab (RAM)	IGHV3-21*01	IGKV1-12*01	90.7	66.0	900.0	9.4	8.7	0.00	0.0
Ustekinumab (UST)	IGHV5-51*01	IGKV1D-16*01	152.7	69.5	1000.0	8.8	8.6	0.15	0.4

Summary

- Rapid Discovery Workflow:** Invenra has developed an accelerated process for discovering novel bispecific antibodies, delivering comprehensive lead panels within just four months.
- Optimized Phage Libraries:** The process begins with phage libraries engineered to maximize diversity while maintaining exceptional humanness and developability.
- Seamless Lead Development:** Identified lead binders are rapidly expressed, characterized, and formatted into the B-Body Platform for further evaluation.
- Streamlined Evaluation:** The B-Body Platform enables high-throughput assessment of all binder permutations, achieving high antibody yields and single-pass purity through an efficient workflow.
- Scalable Production:** The platform supports high-yield production in both transient and stable engineered cell lines, with observed yields reaching up to 11 g/L.

References

- Protein Engineering, Design and Selection, Volume 32, Issue 7, July 2019, Pages 347–354 Bioinformatics, Volume 37, Issue 22, November 2021, Pages 4041–4047
- Jain T, Sun T, Durand S, *et al.* Biophysical properties of the clinical-stage antibody landscape. *Proc Natl Acad Sci U S A.* 2017;114(5):944-949.

Example Transient Yield of a B-Body Lead

Expression Titer	442.8 mg/L
Purification Resin	MabSelect Sure
Molecular Weight	146270 Da
Isoelectric Point	7.8

Biophysical Properties	PDI	Z-Ave (nm)	Tm (°C)	Tagg (°C)
	0.410	71.6	68.7	66.2




Figure 4. High-Yield of a B-Body Bispecific Lead

A B-Body Bispecific lead was transiently expressed in 100 mL of CHO cells and purified using a single-step Protein A chromatography process. The resulting antibody achieved a purity level exceeding 93%, as confirmed by both SDS-PAGE and HPLC analysis, demonstrating the platform's efficiency in generating high-quality bispecific candidates.

